

BET Application Form

Sample preprocessing : Put in Oven (105°C) over 12 hr to remove water.

Sample weight prepared exceed 0.2g.

Applicant : _____ Date : ____ / ____ / ____ Student ID no : _____

Contact Phone: _____ Advisor Sign : _____

Way of operating

Self-using (GIEE student only, NT. 50/hr)

Entrusting

	GIEE	Other departments NTU	party(ies) not affiliated with NTU
BET	<input type="checkbox"/> NT 1,200 /sample	<input type="checkbox"/> NT 1,700 元/ sample	<input type="checkbox"/> NT 2,200 元/ sample
Mesopore	<input type="checkbox"/> NT 2,200 元/sample	<input type="checkbox"/> NT 3,000 元/ sample	<input type="checkbox"/> NT 3,500 元/ sample
Micropore	<input type="checkbox"/> NT 3,200 元/ sample	<input type="checkbox"/> NT 4,000 元/ sample	<input type="checkbox"/> 4 NT,500 元/ sample

Description of sample :

	Sample 1	Sample 2	Sample3	Sample 4
Sample source				
Manufacturing process				
Estimate BET				
Degas Temperature				
Degas Time				